



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Hospital Bedside Locker DX

**Product Code :**  
SAMCOMEDICAL10021A002

**Description :**

Hospital Bedside Locker DX

**Technical Specification :**

Complete locker box made of 22G CRC machine pressed sheet closed from three sides and a suitable drawer for extra storage. Stainless steel top of 22G in matt finish. Complete locker mounted on two PVC Stumps in front and two 5 cms. dia castors at the back for easy moving. Pre-treated and epoxy powder coated finish. Size(in) : 16 L x 16 W x 32 H Hospital Bedside Locker DX, Bedside Lockers Exporters, Bedside Lockers Manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),