



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Surgical Tape, Hypo-Allergenic Microporous

**Product Code :**  
SAMCOMEDICAL500001A127



**Description :**

Surgical Tape, Hypo-Allergenic Microporous Surgical Tape is backed by a Thin highly permeable micro porous non woven fabric. Uniformly coated with hypo allergenic adhesive. Safe to use on sensitive part of the body. Suitable for infants and people with allergies. Peel off easily without pulling out hair or causing pain. High resistance to water and vapour and therefore will remain in position when immersed in water



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),