



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Disposable Reinforced Silicone Laryngeal Mask Airway

**Product Code :**  
SAMCOMEDICAL500001A134



**Description :**

Disposable Reinforced Silicone Laryngeal Mask Airway 100% Medical Grade Silicone. A Stainless Steel spring is inserted in the tube to reduce the risk of kinking



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),