



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Reusable Reinforced Silicone Laryngeal Mask Airway

Product Code :
SAMCOMEDICAL500001A135



Description :

Reusable Reinforced Silicone Laryngeal Mask Airway 100% Medical Grade Silicone. A Stainless Steel spring is inserted in the tube to reduce the risk of kinking. Can be autoclaved repeatedly at 134°C. Reusable up to 40 times



Samco Medical

Website: www.samcomedical.com,