



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Blood Bag Systems

Product Code :
SAMCOMEDICAL500001A146



Description :

Blood Bag Systems Siliconised ultra thin walled, 16G, sharp Japanese needle assures smooth and atraumatic venipuncture. Soft, frosted and flexible kink resistant tubing with unique identification number at intervals of 10 cm Tamper proof, safe and easy to open port covers to prevent contamination Manufactured from high quality medical grade PVC sheet Individually packed in special multi layer pouch and multiples packed in moisture barrier aluminum pouch Convenient hanger slits and holes allow easy suspension of the bags during collection & storage Anticoagulant CPD solution preserves whole blood up to 21 days Anticoagulant CPDA solution preserves whole blood up to 35 days SAGM solution preserves RBC upto 42 days



Samco Medical

Website: www.samcomedical.com,