



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Blood Administration Set with Airvent

**Product Code :**  
SAMCOMEDICAL500001A147



**Description :**

Blood Administration Set with Airvent For transfusion of blood or blood components Transparent & flexible drip chamber Approximately 20 drops/ml Soft and kink resistant PVC tubing With 200 micron blood filter



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),