



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Samcomedical Semi Fowler Bed Manual

**Product Code :**  
SAMCOMEDICAL1001A028

**Description :**

Samcomedical Semi Fowler Bed Manual

**Technical Specification :**

Samcomedical Semi Fowler Bed Manual with mechanically operated back rest by super-smooth crank mechanism. Size: L 2140mm

W 940mm

H 500mm. Removable & interchangeable high quality beautiful ABS engineering plastic head panel and foot panel. Collapsible side railings. (Top bar Aluminium, Vertical bars SS, Bottom pipe epoxy) Epoxy coated mild steel frame work and 2 section perforated top. 125mm dia noiseless castors (all with individual brakes) Provision for I.V. Rod on both sides of the bed. Freight saving knock-down construction



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),