



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Standard Semi-Fowler Bed

Product Code :
SAMCOMEDICAL1001A031

Description :

Standard Semi-Fowler Bed

Technical Specification :

Standard Semi-Fowler Bed Standard semi-fowler bed with mechanically operated back rest by crank mechanism. Size: 2030 $\frac{1}{2}$ 900 $\frac{1}{2}$ 530 mm. Epoxy coated mild steel frame work and 2 section top. Tubular head & foot bows made of stainless steel. Mounted on protective stumps. Provision for I.V. Rod. Freight saving knock-down construction



Samco Medical

Website: www.samcomedical.com,