



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Standard Semi-Fowler Beds

**Product Code :**  
SAMCOMEDICAL1001A032

**Description :**

Standard Semi-Fowler Bed

**Technical Specification :**

Standard Semi-Fowler Bed Standard semi-fowler bed with mechanically operated back rest by crank mechanism. Size: 2030  $\frac{1}{2}$  900  $\frac{1}{2}$  500 mm. Epoxy coated mild steel frame work and 2 section top. Mounted on protective stumps. Provision for I.V. Rod. Freight saving knock-down construction



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),