



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Semi-Fowler Bed (Deluxe)

Product Code :
SAMCOMEDICAL1001A033

Description :

Semi-Fowler Bed (Deluxe)

Technical Specification :

Semi-Fowler Bed (Deluxe) Size: 206L x 90W x 60H cms. Frame work made of Rectangular M.S. Tube 2 Section Top made of perforated M.S. Sheet. Epoxy Powder Coated Finish. Back Rest Section, manoeuvred by Screw Handle from Foot End Detachable S.S. Head and Foot Bows fitted with laminated panels Provision for four saline locations. Legs fitted with PVC Stumps Knockdown Construction



Samco Medical

Website: www.samcomedical.com,