



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Semi-Fowler Bed (General)

**Product Code :**  
SAMCOMEDICAL1001A035

**Description :**

Semi-Fowler Bed (General)

**Technical Specification :**

Semi-Fowler Bed (General) Size: 206L x 90W x 60H cms. Frame work made of Rectangular M.S. Tube. Epoxy Powder Coated Finish. 2 Section Top made of perforated CRCA Sheet with mattress guards. Back Rest Section, manoeuvred by Screw Handle from Foot End Tubular head and foot bows of unequal height. Provision for four saline locations. Legs fitted with PVC Stumps Knockdown Construction



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),