



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Hospital Bed (two Section), DELUXE (With ABS Head & Foot panels)

Product Code :
SAMCOMEDICAL1001A036

Description :

Hospital Bed (two Section), DELUXE (With ABS Head & Foot panels)

Technical Specification :

Hospital Bed (two Section), DELUXE (With ABS Head & Foot panels) Back rest section maneuvered smoothly by screw mechanism (hand crank) from foot end. Overall size : 2030L x 900W x 600H mm. Frame work of mild steel rectangular CRC tubes. Mounted on protective stumps. Two section perforated CRC sheet top. Head & Foot panels made of ABS high quality. Provision for saline / IV rod at four locations. Pre-treated and epoxy powder coated finish. Knockdown construction



Samco Medical

Website: www.samcomedical.com,