



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :

Hospital Bed Two Section, General (With built in back rest and rigid mesh mattress platform)

Product Code :

SAMCOMEDICAL1001A037

Description :

Hospital Bed Two Section, General (With built in back rest and rigid mesh mattress platform)

Technical Specification :

Hospital Bed Two Section, General (With built in back rest and rigid mesh mattress platform) Back rest adjustable on lever mechanism. Overall size : 203L x 90W x 60H cms. Frame work of mild steel rectangular CRC tubes. Mounted on PVC stumps. Head and leg bow of M.S. tubular pipes, epoxy coated with vertical supports on either side. Provision for saline / IV rod at four locations. Pre-treated epoxy powder coated finish. Knockdown construction



Samco Medical

Website: www.samcomedical.com,