



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Semi-Fowler Bed, Manual (with Provision for Bed Pan)

Product Code :
SAMCOMEDICAL1001A038

Description :

Semi-Fowler Bed, Manual (with Provision for Bed Pan)

Technical Specification :

Semi-Fowler Bed, Manual (with Provision for Bed Pan) Semi-Fowler Bed with mechanically operated back rest by crank mechanism. Size: 2000 ĩ½ 900 ĩ½ 500 mm. Provision for Bed pan. Epoxy coated mild steel frame work and 2 section top. Mounted on protective stumps. Provision for I.V. Rod on both sides of the bed. Freight saving knock-down construction



Samco Medical

Website: www.samcomedical.com,