



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Samcomedical Fowler Bed Manual

Product Code :
SAMCOMEDICAL1001A039

Description :

Samcomedical Fowler Bed Manual

Technical Specification :

Samcomedical Fowler Bed Manual with mechanically operated back rest and knee rest by super-smooth crank mechanism. Size: L 2140mm

W 940mm

H 500mm. Two easy lifting guardrails (1 on each side), which are safe and reliable, and can be fixed upward and downward. Removable & interchangeable high quality beautiful ABS engineering plastic head panel and foot panel. Epoxy coated mild steel frame work and 4 section perforated top. 125mm dia noiseless castors with Central braking system. Provision for I.V. Rod



Samco Medical

Website: www.samcomedical.com,