



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Samcomedical Fowler Bed Manual With Circulated

Product Code :
SAMCOMEDICAL1001A041

Description :

Samcomedical Fowler Bed Manual With Circulated Central Braking System

Technical Specification :

Samcomedical Fowler Bed Manual With Circulated Central Braking System with mechanically operated back rest and knee rest by super-smooth crank mechanism. Size: L 2140mm

W 940mm

H 530mm. Four easy lifting guardrails (2 on each side), which are safe and reliable, and can be fixed upward and downward. Removable & interchangeable high quality beautiful ABS engineering plastic head panel and foot panel. D type design. Epoxy coated mild steel frame work and 4 section perforated top. Circulated Central braking system locks / unlocks all swivel castor simultaneously. Provision for I.V. Rod



Samco Medical

Website: www.samcomedical.com,