



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Samcomedical Fowler Bed Manual Two Function

Product Code :
SAMCOMEDICAL1001A044

Description :

Samcomedical Fowler Bed Manual Two Function

Technical Specification :

Samcomedical Fowler Bed Manual Two Function with mechanically operated back rest and knee rest by super-smooth crank mechanism. Size: L 2140mm

W 940mm

H 500mm. Removable & interchangeable high quality beautiful ABS engineering plastic head panel and foot panel. Collapsible side railings. (Top bar Aluminium, Vertical bars SS, Bottom pipe epoxy) Epoxy coated mild steel frame work and 4 section perforated top. 125mm dia noiseless castors (all with individual brakes) Provision for I.V. Rod



Samco Medical

Website: www.samcomedical.com,