



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Samcomedical Fowler Bed Manual 2 Function with Food Tray

**Product Code :**  
SAMCOMEDICAL1001A047

**Description :**

Samcomedical Fowler Bed Manual 2 Function with Food Tray

**Technical Specification :**

Samcomedical Fowler Bed Manual 2 Function with Food Tray with mechanically operated back rest and knee rest by super-smooth crank mechanism. Size: L 2170mm

W 940mm

H 500mm. Removable & interchangeable high quality beautiful ABS engineering plastic head panel and foot panel. Collapsible side railings. (Top bar Aluminium, Vertical bars SS, Bottom pipe epoxy) Equipped with Food Tray. Epoxy coated mild steel frame work and 4 section perforated top. 125mm dia noiseless castors (all with individual brakes) Provision for I.V. Rod on both sides of the bed. Freight saving knock-down construction



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),