



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Hospital Bed(four Section), deluxe(With ABS Head & Foot panels)

**Product Code :**  
SAMCOMEDICAL1001A051

**Description :**

Hospital Bed(four Section), deluxe(With ABS Head & Foot panels)

**Technical Specification :**

Hospital Bed(four Section), deluxe(With ABS Head & Foot panels) Back rest and knee rest are maneuvered smoothly by separate screw mechanism (hand crank) from the foot end. Overall size : 2030L x 900W x 600H mm. Head & Foot Panels made of ABS high quality. Frame work of mild steel rectangular CRC tubes. Mounted on 10 cms dia heavy duty castors (2 with brakes). Four sections perforated CRC sheet top. Provision for saline / IV rod at four locations. Pre-treated and epoxy powder coated finish. Knockdown construction



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),