



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Hospital Foot Stool Double Step

**Product Code :**  
SAMCOMEDICAL10016A002

**Description :**

Hospital Foot Stool Double Step

**Technical Specification :**

Hospital Foot Stool Double Step: Size 51L x 30W x 23/45H cms, ERW steel tubular frame work mounted on PVC shoes. Both step fixed with rubber mat surrounded by aluminium angle frame. Pre-treated and epoxy powder coated.



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),