



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Hospital Foot Stool Double Step

Product Code :
SAMCOMEDICAL10016A002

Description :

Hospital Foot Stool Double Step

Technical Specification :

Hospital Foot Stool Double Step: Size 51L x 30W x 23/45H cms, ERW steel tubular frame work mounted on PVC shoes. Both step fixed with rubber mat surrounded by aluminium angle frame. Pre-treated and epoxy powder coated.



Samco Medical

Website: www.samcomedical.com,