



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Samcomedical Standard Fowler Provision Bed

Product Code :
SAMCOMEDICAL1001A054

Description :

Samcomedical Standard Fowler Provision Bed

Technical Specification :

Standard Fowler Bed Standard fowler bed with mechanically operated back rest and knee rest by crank mechanism. Size: 2030 900 500 mm. Epoxy coated mild steel frame work & 4 section top.. Mounted on protective stumps. Provision for I.V. Rod. Freight saving knock-down construction



Samco Medical

Website: www.samcomedical.com,