



Email : sales@samcomedical.com  
Phone: 1-300-400-8211

**Product Name :**  
Fowler Bed (General)

**Product Code :**  
SAMCOMEDICAL1001A056

**Description :**

Fowler Bed (General)

**Technical Specification :**

Fowler Bed (General) Back rest and knee rest sections maneuvered smoothly by separate screw mechanism (hand crank) from foot end. Overall size : 2030L x 900W x 600H mm. Frame work of mild steel rectangular CRC tubes. Tubular head and foot bows of unequal height with vertical tube supports mounted on 10 cms heavy duty castors (2 with brakes). Four sections perforated CRC sheet top. Provision for saline / IV rod at four locations. Pre-treated and epoxy powder coated finish. Knockdown construction



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),