



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Fowler Bed (Deluxe)

Product Code :
SAMCOMEDICAL1001A057

Description :

Fowler Bed (Deluxe)

Technical Specification :

Fowler Bed (Deluxe) Back rest and knee rest are maneuvered smoothly by separate screw mechanism (hand crank) from the foot end. Overall size : 2030L x 900W x 600H mm. Frame work of mild steel rectangular CRC tubes. "H" type legs on 10 cms dia heavy duty castors (2 with brakes). Four sections perforated CRC sheet top. Stainless steel head & foot bow panels. Provision for saline / IV rod at four locations. Pre-treated and epoxy powder coated finish. Knockdown construction



Samco Medical

Website: www.samcomedical.com,