



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Fowler Beds (Deluxe)

Product Code :
SAMCOMEDICAL1001A058

Description :

Fowler Beds (Deluxe)

Technical Specification :

Fowler Bed (Deluxe) Size: 206L x 90W x 60H cms. Frame work made of Rectangular M.S. Tube 4 Section Top made of perforated M.S. Sheet. Epoxy Powder Coated Finish. Back Rest, Knee Rest Positions manoeuvred by Separate screw from Foot End. Detachable S.S. Head and Foot Bows fitted with laminated panels Telescopic IV Rod with four Locations. Bed Mounted on 100mm Dia Castor- two with brakes. Knockdown Construction



Samco Medical

Website: www.samcomedical.com,