



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :

Samcomedical Fowler Bed Manual 2 Function Bed Pan

Product Code :

SAMCOMEDICAL1001A060

Description :

Samcomedical Fowler Bed Manual 2 Function with Provision for Bed Pan

Technical Specification :

Fowler Bed, Manual, Two Function (with Provision for Bed Pan) Fowler Bed with mechanically operated back rest 0 to 75 and knee rest 0 to 40 by crank mechanism. Size: 2000 980 500 mm. Provision for Bed pan. Collapsible side railings. ABS Head & Foot panels. Epoxy coated mild steel frame work and 4 section top. Mounted on protective stumps. Provision for I.V. Rod on both sides of the bed. Freight saving knock-down construction



Samco Medical

Website: www.samcomedical.com,