



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Samcomedica Fowler Bed Manual with Bed Pan

**Product Code :**  
SAMCOMEDICAL1001A062

**Description :**

Samcomedical Fowler Bed Manual with Provision for Bed Pan

**Technical Specification :**

Fowler Bed, Manual (with Provision for Bed Pan) Fowler Bed with mechanically operated back rest 0 to 75 and knee rest 0 to 40 by crank mechanism. Size: 2000 900 500 mm. Provision for Bed pan. With side railings. Epoxy coated mild steel frame work and 4 section top. Mounted on protective stumps. Provision for I.V. Rod on both sides of the bed. Freight saving knock-down construction



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),