



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :

Samcomedical Standard Plain Steel Hospital Bed

Product Code :

SAMCOMEDICAL1001A065

Description :

Samcomedical Standard Plain Steel Hospital Bed

Technical Specification :

Standard Plain Hospital Bed Size: 2030 900 500 mm. Epoxy coated mild steel frame work,with mesh wire platform. Provision for I.V. Rod. Freight saving knock-down construction with castor



Samco Medical

Website: www.samcomedical.com,