



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :

Samcomedical Hospital Bed Plain

Product Code :

SAMCOMEDICAL1001A066

Description :

Samcomedical Hospital Bed Plain General With rigid mesh mattress platform

Technical Specification :

Hospital Bed Plain, General (With rigid mesh mattress platform) Overall size : 2030L x 900W x 600H mm Frame work of mild steel rectangular CRC tubes. Mounted on heavy duty castors (2 with brakes). Head and leg bow of M.S. tubular pipe vertical supports on either side. Provision for saline / IV rod at four locations. Pre-treated epoxy powder coated finish. Knockdown construction



Samco Medical

Website: www.samcomedical.com,