



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Attendant Bed

Product Code :
SAMCOMEDICAL1001A067

Description :

Attendant Bed

Technical Specification :

Attendant Bed Size: 180L x 60W x 46H cms. Frame work made of rectangular / square M.S. Tube. Epoxy Powder Coated Finish 5 Cm Foam Padded top covered with rexine. Legs fitted with PVC stumps. Knockdown Construction. * CRCA : Cold Rolled Closed Annealed * M.S. : Mild Steel



Samco Medical

Website: www.samcomedical.com,