



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Orthopaedic Bed, Fowler

Product Code :
SAMCOMEDICAL1001A068

Description :

Orthopaedic Bed, Fowler

Technical Specification :

Orthopaedic Bed, Fowler Orthopaedic bed with mechanically operated back rest & knee rest all by Super Smooth crank mechanism. Size: 2150 $\frac{1}{2}$ 940 $\frac{1}{2}$ 500 mm. Epoxy Coated Steel frame work and 4 section perforated top. Stainless Steel Orthopaedic attachment (without weight chain and weight hangers). Head panel and foot panel equipped with safety lock and roller bumpers. Mounted on protective stumps. Provision for I.V. Rod. Freight saving knock-down construction



Samco Medical

Website: www.samcomedical.com,