



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Orthopaedic Bed, with ABS Panels

Product Code :
SAMCOMEDICAL1001A069

Description :

Orthopaedic Bed, with ABS Panels

Technical Specification :

Orthopaedic Bed, with ABS Panels Orthopaedic bed with mechanically operated back rest & knee rest by super-smooth crank mechanism. Size: 2150 $\frac{1}{2}$ 940 $\frac{1}{2}$ 500 mm. Epoxy coated mild steel frame work and 4 section perforated top. Removable high quality ABS head & foot panels. Mounted on protective stumps. Bumper Protection at four corners. Provision for I.V. Rod. Freight saving knock-down construction



Samco Medical

Website: www.samcomedical.com,