



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Orthopadic Bed (With SS Frame & ABS Panel)

Product Code :
SAMCOMEDICAL1001A071

Description :

Orthopadic Bed (With SS Frame & ABS Panel)

Technical Specification :

Orthopadic Bed (With SS Frame & ABS Panel) Orthopedic traction bed with cranks operated back rest tilting 0-70°, Knee rest tilting 0-40°, Height Adjustment (730mm to 530mm). Left and right leg part can be raised up separately, convenient for treatment and use. Size: L 2235mm x W 960mm x H 530-730mm (Adjustable Height). Epoxy coated Mild steel frame work and 4 section perforated top. Removable High Quality ABS engineering Plastic head & foot panel. Size: 1910mm x 830mm Collapsible aluminum side rails and stainless steel frame with pulley. 125mm dia noiseless castors, 2 with brakes



Samco Medical

Website: www.samcomedical.com,