



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Orthopaedic Bed, S.S. Multi Function

Product Code :
SAMCOMEDICAL1001A072

Description :

Orthopaedic Bed, S.S. Multi Function

Technical Specification :

Orthopaedic Bed, S.S. Multi Function Orthopaedic bed with mechanically operated back rest 0 to 70°½ & knee rest 0 to 35°½ and trendelenburg (TR) all by crank mechanism. Knee-rest portion sectioned in two parts for patient's convenience. Size: 2030 x 900 x 530 mm. Stainless steel frame work and 5 section top. Balkan Frame with monkey bars and monkey chains. Pulley attachment (without weight chain and weight hangers). Four swivel castors (all with individual brakes). Provision for I.V. Rod. Freight saving knock-down construction



Samco Medical

Website: www.samcomedical.com,