



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Ambulance Stretcher

Product Code :
SAMCOMEDICAL10025A001

Description :

Ambulance Stretcher

Technical Specification :

Light weight, strong structure, lovely design, good arts and crafts, easy antiseptis, small size, convenient to use. Utilizes high precision mechanical transmission. Can be used automatically folded through the left and right hand control handle. Made of high quality Aluminium alloy tubes. Applicable for various kind of hospitals, emergency centers, ambulances and battle fields. Dimensions: A - High position : 1,900 x 560 x 900mm B - Low position: 1,900 x 560 x 280mm Load weight: > = 120kg Product weight: 40 kg Diameter of the castors: 125mm
Ambulance Stretcher, Hospital Stretchers Exporters, Hospital Stretchers Manufacturers



Samco Medical

Website: www.samcomedical.com,