



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Infant Bed / Child Cot / Baby Bassinet

Product Code :
SAMCOMEDICAL10029A002

Description :

Infant Bed / Child Cot / Baby Bassinet

Technical Specification :

Infant Bed / Child Cot / Baby Bassinet Size: L 105mm i 1/2 W 45mm i 1/2 H 1045mm i 1/2 855mm (Height Adjustable by hydraulic system). In build touch keyboard. Weight scale to Measure baby weight . Tilt angle 15i 1/2, easy for mother and baby to communicate and feed. Hand rail could be placed flat or at 90i 1/2 angle. Can be used as a mosquito net support. Additional Features: SD memory card to play music, can store music in MP3 format Atmosphere temperature measuring



Samco Medical

Website: www.samcomedical.com,