



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Swaddling Table

**Product Code :**  
SAMCOMEDICAL10029A003

**Description :**

Swaddling Table

**Technical Specification :**

Swaddling Table Intended for babies swaddle. Upper part covered by foam-rubber and padded by leatherette. Height: 980mm.  $\pm 10$ mm. Length: 850mm.  $\pm 10$ mm. Width: 700mm.  $\pm 10$ mm. Framework made of mild steel epoxy coated finish. The table is easy for disinfecting



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),