



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Infant Bed / Child Cot, S.S.

Product Code :
SAMCOMEDICAL10029A004

Description :

Infant Bed / Child Cot, S.S.

Technical Specification :

Infant Bed / Child Cot, S.S. Size: 810 ½ 460 ½ 960 mm. Completely made of stainless steel. 1 S.S. shelf at the bottom & S.S. guard rails. 4 swivel castor (2 with brakes). Freight saving knock-down construction



Samco Medical

Website: www.samcomedical.com,