



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :

Infant Bed / Child Cot, Epoxy Coated Steel

Product Code :

SAMCOMEDICAL10029A005

Description :

Infant Bed / Child Cot, Epoxy Coated Steel

Technical Specification :

Infant Bed / Child Cot, Epoxy Coated Steel Size: 810 ½ 460 ½ 960 mm. Completely made of Epoxy Coated Steel. 1 shelf at the bottom & guard rails made of epoxy coated steel. 4 swivel castor (2 with brakes). Freight saving knock-down construction



Samco Medical

Website: www.samcomedical.com,