



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Ambulance Foldaway Stretcher

Product Code :
SAMCOMEDICAL10025A002

Description :

Ambulance Foldaway Stretcher

Technical Specification :

Ambulance Stretcher Foldaway Stretcher with wheels, outriggers and can be folded. It is made of high strength aluminum alloy material and the waterproof leather materials. It has characterizes of light-weighted, small-sized, easy-carried, reliable quality, use-safely and easy for sterilization. Dimensions Product size (L x W x H): 185 x 50x 20cm Folded size (L x W x H): 92.5x 50 x 9cm Packing size (2pcs/carton): 96 x 54 x 16cm N.Weight: 6kg G.W.:7 kg Load bearing: less than 159kg Ambulance Foldaway Stretcher, Hospital Stretchers Exporters, Hospital Stretchers Manufacturers



Samco Medical

Website: www.samcomedical.com,