



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Semi-Fowler Bed For Children, with Side Railings

Product Code :
SAMCOMEDICAL10029A008

Description :

Semi-Fowler Bed For Children, with Side Railings

Technical Specification :

Semi-Fowler Bed For Children, with Side Railings Semi-Fowler Bed with mechanically operated Back rest by crank mechanism. Size: 1880 $\frac{1}{2}$ 880 $\frac{1}{2}$ 580 mm. Epoxy coated mild steel frame work with 2 section top. S.S. side railings with Up & Down movement. With S.S. Shelf at the bottom for keeping blankets, basin etc. Tubular Head & Foot Bows also made of stainless steel. Mounted on four silent swivel castor (all with brakes). Provision for I.V. Rod. Freight saving knock-down construction



Samco Medical

Website: www.samcomedical.com,