



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Samcomedical Semi-Fowler Bed for Children

**Product Code :**  
SAMCOMEDICAL10029A010

**Description :**

Samcomedical Semi-Fowler Bed for Children, with Side Railings

**Technical Specification :**

Semi-Fowler Bed for Children, with Side Railings Semi-Fowler Bed with mechanically operated Back rest by Super-Smooth crank mechanism. Size: 1650 825 630 mm. Epoxy coated mild steel frame work with 2 section top. Dropdown side railings. Mounted on four silent swivel castor (2 with brakes). Provision for I.V. Rod. Freight saving knock-down construction



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),