



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**

Baby Bassinet Deluxe (With poly crib)

**Product Code :**

SAMCOMEDICAL10029A011

**Description :**

Baby Bassinet Deluxe (With poly crib)

**Technical Specification :**

Baby Bassinet Deluxe (With poly crib) Frame work made of ERW tubular pipe of 2.5 cms. dia. Mounted on four 5 cms dia. castors. Removable CRIB of poly glass tiltable to trendelenburg and R-trend positions. Laminated bottom shelf with a spacious cabinet for storage. Provision for monitor stand and drip rod available. Pre-treated and epoxy powder coated



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),