



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Samcomedical Baby Cradle On Stand Swinging Type

**Product Code :**  
SAMCOMEDICAL10029A012

**Description :**

Samcomedical Baby Cradle On Stand Swinging Type

**Technical Specification :**

Baby Cradle On Stand Swinging Type Size : 36 x 15W x 39H inches. Height upto CRIB top 98 cms. Stand made of ERW tubular pipe with hanging hooks mounted on rubber shoes. CRIB made of steel bars/strips. Pre-treated and epoxy powder coated



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),