



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Pediatric Beds (With Side Railings)

**Product Code :**  
SAMCOMEDICAL10029A013

**Description :**

Pediatric Beds (With Side Railings)

**Technical Specification :**

Pediatric Bed (With Side Railings) Frame work of rectangular CRC material. Mounted on four 10 cms. dia castors (2 with brakes). Drop side, safety railing at both side. Pre-treated and epoxy powder coated finish. Knockdown construction



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),