



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Overbed Table (Adjustable With Geared Handle)

**Product Code :**  
SAMCOMEDICAL10014A031



**Description :**

Overbed Table (Adjustable With Geared Handle) Table base and frame work made of strong rectangular / square steel tubes mounted on four 5 cms dia. castors. Overall size of laminated top : 76L x 40W cms. adjustable by gear handle. Pre-treated and epoxy powder coated. Knockdown construction



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),