



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Ambulance Foldaway Stretcher DX

Product Code :
SAMCOMEDICAL10025A003

Description :

Ambulance Foldaway Stretcher DX

Technical Specification :

Ambulance Foldaway Stretcher DX with wheels, outriggers and can be folded. It is made of high strength aluminum alloy and the water-tightness leather material. The back of this stretcher can be lifted, and the maximum angle is 60°. It is characterized by its being fashionable structure, light-weighted, easy-carried, use-safely, reliable quality and easy for sterilization. Dimensions Product size (L x W x H): 185 x 50 x 24cm Folding size (L x W x H): 92.5 x 50 x 11cm Packing size (1pcs): 96 x 54 x 12cm Ambulance Foldaway Stretcher DX, Hospital Stretchers Exporters, Hospital Stretchers Manufacturers



Samco Medical

Website: www.samcomedical.com,