



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Foldaway Stretcher Aluminum

Product Code :
SAMCOMEDICAL10025A004

Description :

Foldaway Stretcher Aluminum

Technical Specification :

Foldaway Stretcher Aluminum : This folding stretcher is made of high strength aluminum alloy material. It has characterizes of small-sized, easy-carried, reliable quality, use-safely and easy for sterilization, etc. It is commonly used in hospital, gymnasium, ambulance, and army hospital to carry patients and wounded person. Furthermore, it is much more used for large acreage fractured patients, and the hole of the back and front of this stretcher can be used with head immobilize at the same time. Dimensions Product size (L x W x H):185 x 50 x 3cm Folding size (L x W x H):92.5 x 50 x 7cm Packing size (2pcs/carton): 96 x 54 x16cm N.W.: 6.5kg G.W.: 7.6kg Load bearing: less than 159kg Foldaway Stretcher Aluminum, Hospital Stretchers Exporters, Hospital Stretchers Manufacturers



Samco Medical

Website: www.samcomedical.com,