



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Foldaway Stretcher Army

**Product Code :**  
SAMCOMEDICAL10025A005

**Description :**

Foldaway Stretcher Army

**Technical Specification :**

Foldaway Stretcher Army: This stretcher is made of high strength aluminum alloy material and the oxford leather materials. It has characterizes of light-weighted, small-sized, easy-carried, reliable quality, use-safely and easy for sterilization, etc. It can be folded by transverse. It is mainly used for battle field and outdoor carrying patients and wounded person. Furthermore, the hole of the back and front of this stretcher can be used with head immobilize at the same time. Dimensions: Product size (L x W x H): 185 x 48 x 3cm Folded size (L x W x H): 46 x 48 x 14cm Packing size (2pcs/carton): 50x 48 x 15cm,br> 4) N.W.: 4.5kg G.W.: 6.2kg Load bearing: less than 159kg Foldaway Stretcher Army, Hospital Stretchers Exporters, Hospital Stretchers Manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),