



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Folding Stretcher

Product Code :
SAMCOMEDICAL10025A006

Description :

Folding Stretcher

Technical Specification :

Folding Stretcher with outrigger is made of high strength aluminum alloy material and the waterproof leather materials. This stretcher with outrigger can be folded in length and width, easy for transportation and storage. The orange fabric is fire proof, water proof, anti-cracking and easy to clean. It is characterized by its being light weight, portable, reliable quality, easy to carry out and safe to use, etc. It is mainly used for the hospitals, emergency centers and battle field to carry patients and wounded person. Dimensions: Product size (L x W x H): 208 x 54.5 x 13cm Folding size (L x W x H): 104 x 9 x 17cm Packing size (2pcs): 108 x 18 x 18.5cm N.W.: 4.7kg G.W.: 5.3kg Load bearing: less than 159kg Folding Stretcher, Hospital Stretchers Exporters, Hospital Stretchers Manufacturers



Samco Medical

Website: www.samcomedical.com,